

Product Datasheet

Mouse ATP1b4 ELISA Kit (orb553022)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	27.8pg/mL
Expiration Date	Please enquire.