

Product Datasheet

Mouse BLC1 ELISA Kit (Ready to Use) (orb553031)

Description	Mouse B-Lymphocyte Chemoattractant 1 (BLC1) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O55038
Sensitivity	5.9pg/mL
Entrez	55985
Expiration Date	Please enquire.