

Product Datasheet

Mouse bTG ELISA Kit (Ready to Use) (orb553037)

Description	Mouse Beta-Thromboglobulin (bTG) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P62325
Sensitivity	34pg/mL
Entrez	12226
Expiration Date	Please enquire.