

Product Datasheet

Mouse ROS1 ELISA Kit (Ready to Use) (orb553061)

Description	Mouse C-Ros Oncogene 1, Receptor Tyrosine Kinase (ROS1) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	33pg/mL
Expiration Date	Please enquire.