

Product Datasheet

Mouse CNP ELISA Kit (Ready to Use) (orb553062)

Description	Mouse C-Type Natriuretic Peptide (CNP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P16330
Sensitivity	12.6pg/mL
Entrez	12799
Expiration Date	Please enquire.