

Product Datasheet

Mouse CACNa1B ELISA Kit (Ready to Use) (orb553074)

Description	Mouse Calcium Channel, Voltage Dependent, N-Type, Alpha 1B Subunit (CACNa1B) Sandwich ELISA KitReady...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	055017
Sensitivity	0.053ng/mL
Entrez	12287
Expiration Date	Please enquire.