

Product Datasheet

Mouse PCDHA6 ELISA Kit (orb553082)

Description	Mouse Cannabinoid Receptor 2, Macrophage (PCDHA6) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.12-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	30pg/mL
Expiration Date	Please enquire.