

Product Datasheet

Mouse CAT2 ELISA Kit (Ready to Use) (orb553121)

Description	Mouse Cationic Amino Acid Transporter 2 (CAT2) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.124ng/mL
Expiration Date	Please enquire.