

Product Datasheet

Mouse CEBPb ELISA Kit (Ready to Use) (orb553123)

Description	Mouse CCAAT/Enhancer Binding Protein Beta (CEBPb) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28033
Sensitivity	33pg/mL
Entrez	12608
Expiration Date	Please enquire.