

Product Datasheet

Mouse CCR10 ELISA Kit (Ready to Use) (orb553136)

Description	Mouse Chemokine C-C-Motif Receptor 10 (CCR10) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JL21
Sensitivity	0.051ng/mL
Entrez	12777
Expiration Date	Please enquire.