

Product Datasheet

Mouse CCR7 ELISA Kit (Ready to Use) (orb553138)

Description	Mouse Chemokine C-C-Motif Receptor 7 (CCR7) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P47774
Sensitivity	0.054ng/mL
Entrez	12775
Expiration Date	Please enquire.