

Product Datasheet

Mouse CXCR4 ELISA Kit (Ready to Use) (orb553141)

Description	Mouse Chemokine C-X-C-Motif Receptor 4 (CXCR4) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P70658
Sensitivity	27.2pg/mL
Entrez	12767
Expiration Date	Please enquire.