

## Product Datasheet

### Mouse CHAF1B ELISA Kit (Ready to Use) (orb553153)

|                          |                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Chromatin Assembly Factor 1, Subunit B (CHAF1B) Sandwich ELISA KitReady to Use: Pre-diluted de...           |
| <b>Reactivity</b>        | Mouse                                                                                                             |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                     |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                      |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                               |
| <b>Note</b>              | For research use only                                                                                             |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids. |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.                                                     |
| <b>Uniprot ID</b>        | <b>Q9D0N7</b>                                                                                                     |
| <b>Sensitivity</b>       | 0.052ng/mL                                                                                                        |
| <b>Entrez</b>            | <b>110749</b>                                                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                                                   |