

Product Datasheet

Mouse CBX3 ELISA Kit (Ready to Use) (orb553154)

Description	Mouse Chromobox Homolog 3 (CBX3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P23198
Sensitivity	0.046ng/mL
Entrez	12417
Expiration Date	Please enquire.