

Product Datasheet

Mouse CD34 ELISA Kit (Ready to Use) (orb553169)

Description	Mouse Cluster Of Differentiation 34 (CD34) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q64314
Sensitivity	0.061ng/mL
Entrez	12490
Expiration Date	Please enquire.