

Product Datasheet

Mouse F12 ELISA Kit (Ready to Use) (orb553180)

Description	Mouse Coagulation Factor XII (F12) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityplasma.
Sample Types	plasma.
Uniprot ID	Q80YC5
Sensitivity	0.251ng/mL
Entrez	58992
Expiration Date	Please enquire.