

Product Datasheet

Mouse CART ELISA Kit (Ready to Use) [Out of stock] (orb553182)

Description	Mouse Cocaine And Amphetamine Regulated Transcript (CART) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P56388
Sensitivity	0.055ng/mL
Entrez	27220
Expiration Date	Please enquire.