

Product Datasheet

Mouse Collagen I ELISA Kit (Ready to Use) (orb553184)

Description	Mouse Collagen Type I Alpha 1 (Collagen I) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P11087
Sensitivity	0.137ng/mL
Entrez	12842
Expiration Date	Please enquire.