

Product Datasheet

Mouse CFD ELISA Kit (Ready to Use) (orb553211)

Description	Mouse Complement Factor D (CFD) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Mouse
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P03953
Sensitivity	15.55pg/mL
Entrez	11537
Expiration Date	Please enquire.