

Product Datasheet

Mouse CDK8 ELISA Kit (Ready to Use) (orb553227)

Description	Mouse Cyclin Dependent Kinase 8 (CDK8) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8R3L8
Sensitivity	0.053ng/mL
Entrez	264064
Expiration Date	Please enquire.