

Product Datasheet

Mouse CYP11B1 ELISA Kit (Ready to Use) (orb553241)

Description	Mouse Cytochrome P450 11B1 (CYP11B1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q3TG86
Sensitivity	0.064ng/mL
Entrez	110115
Expiration Date	Please enquire.