

Product Datasheet

Mouse CYP7A1 ELISA Kit (Ready to Use) (orb553244)

Description	Mouse Cytochrome P450 7A1 (CYP7A1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q64505
Sensitivity	0.139ng/mL
Entrez	13122
Expiration Date	Please enquire.