

Product Datasheet

Mouse DLAT ELISA Kit (Ready to Use) (orb553271)

Description	Mouse Dihydrolipoyl Transacetylase (DLAT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8BMF4
Sensitivity	0.048ng/mL
Entrez	235339
Expiration Date	Please enquire.