

Product Datasheet

Mouse DHODH ELISA Kit (Ready to Use) (orb553272)

Description	Mouse Dihydroorotate Dehydrogenase (DHODH) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	035435
Sensitivity	0.056ng/mL
Entrez	56749
Expiration Date	Please enquire.