

Product Datasheet

Mouse DPP9 ELISA Kit (Ready to Use) (orb553275)

Description	Mouse Dipeptidyl Peptidase 9 (DPP9) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8BVG4
Sensitivity	25pg/mL
Entrez	224897
Expiration Date	Please enquire.