

Product Datasheet

Mouse RAD50 ELISA Kit (Ready to Use) (orb553285)

Description	Mouse DNA Repair Protein RAD50 (RAD50) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P70388
Sensitivity	0.113ng/mL
Entrez	19360
Expiration Date	Please enquire.