

Product Datasheet

Mouse DCLK1 ELISA Kit (Ready to Use) (orb553287)

Description	Mouse Doublecortin Like Kinase 1 (DCLK1) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JLM8
Sensitivity	0.059ng/mL
Entrez	13175
Expiration Date	Please enquire.