

Product Datasheet

Mouse EGLN1 ELISA Kit (Ready to Use) (orb553299)

Description	Mouse Egl Nine Homolog 1 (EGLN1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q91YE3
Sensitivity	0.108ng/mL
Entrez	112405
Expiration Date	Please enquire.