

Product Datasheet

Mouse FABP5 ELISA Kit (Ready to Use) (orb553341)

Description	Mouse Fatty Acid Binding Protein 5, Epidermal (FABP5) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, breast milk or other biological fluids.
Sample Types	tissue homogenates, breast milk and other biological fluids.
Uniprot ID	Q05816
Sensitivity	0.054ng/mL
Entrez	16592
Expiration Date	Please enquire.