

Product Datasheet

Mouse FGF3 ELISA Kit (Ready to Use) (orb553365)

Description	Mouse Fibroblast Growth Factor 3 (FGF3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P05524
Sensitivity	5.5pg/mL
Entrez	14174
Expiration Date	Please enquire.