

Product Datasheet

Mouse FGF8 ELISA Kit (Ready to Use) (orb553366)

Description	Mouse Fibroblast Growth Factor 8, Androgen Induced (FGF8) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P37237
Sensitivity	6.6pg/mL
Entrez	14179
Expiration Date	Please enquire.