

Product Datasheet

Mouse FOLR1 ELISA Kit (Ready to Use) (orb553381)

Description	Mouse Folate Receptor 1, Adult (FOLR1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35846
Sensitivity	31pg/mL
Entrez	14275
Expiration Date	Please enquire.