

Product Datasheet

Mouse GAL6 ELISA Kit (Ready to Use) (orb553402)

Description	Mouse Galectin 6 (GAL6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	054891
Sensitivity	0.034ng/mL
Entrez	16857
Expiration Date	Please enquire.