

Product Datasheet

Mouse GCLC ELISA Kit (Ready to Use) (orb553427)

Description	Mouse Glutamate Cysteine Ligase, Catalytic (GCLC) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P97494
Sensitivity	0.118ng/mL
Entrez	14629
Expiration Date	Please enquire.