

Product Datasheet

Mouse GP1Bb ELISA Kit (Ready to Use) (orb553455)

Description	Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) Sandwich ELISA KitReady to Use: Pre-diluted...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityplasma.
Sample Types	plasma.
Uniprot ID	P56400
Sensitivity	14.7pg/mL
Entrez	14724
Expiration Date	Please enquire.