

Product Datasheet

Mouse HSP10 ELISA Kit (Ready to Use) (orb553485)

Description	Mouse Heat Shock 10kDa Protein 1 (HSP10) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q64433
Sensitivity	0.65ng/mL
Entrez	15528
Expiration Date	Please enquire.