

Product Datasheet

Mouse HSP90b1 ELISA Kit (Ready to Use) (orb553493)

Description	Mouse Heat Shock Protein 90kDa Beta 1 (HSP90b1) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08113
Sensitivity	0.064ng/mL
Entrez	22027
Expiration Date	Please enquire.