

## Product Datasheet

### Mouse HSPb6 ELISA Kit (Ready to Use) (orb553494)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Heat Shock Protein Beta 6 (HSPb6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen... |
| <b>Reactivity</b>        | Mouse                                                                                                   |
| <b>Range</b>             | 0.156-10ng/mL                                                                                           |
| <b>Form/Appearance</b>   | Ready to Use                                                                                            |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                     |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.     |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.                                                         |
| <b>Uniprot ID</b>        | <b>Q5EBG6</b>                                                                                           |
| <b>Sensitivity</b>       | 0.067ng/mL                                                                                              |
| <b>Entrez</b>            | <b>243912</b>                                                                                           |
| <b>Expiration Date</b>   | Please enquire.                                                                                         |