

Product Datasheet

Mouse Hepcidin-25 ELISA Kit (Ready to Use) (orb553504)

Description	Mouse Hepcidin (Hepcidin-25) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, urine or other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q9EQ21
Sensitivity	25.4pg/mL
Entrez	84506
Expiration Date	Please enquire.