

## Product Datasheet

### Mouse HAS2 ELISA Kit (Ready to Use) (orb553526)

|                          |                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Hyaluronan Synthase 2 (HAS2)<br>Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...         |
| <b>Reactivity</b>        | Mouse                                                                                                              |
| <b>Range</b>             | 78-5000pg/mL                                                                                                       |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                       |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                                |
| <b>Note</b>              | For research use only                                                                                              |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                                                     |
| <b>Sensitivity</b>       | 26pg/mL                                                                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                                                                    |