

Product Datasheet

Mouse IGFBP7 ELISA Kit (Ready to Use) (orb553557)

Description	Mouse Insulin Like Growth Factor Binding Protein 7 (IGFBP7) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Mouse
Range	6.25-400pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q61581
Sensitivity	2.34pg/mL
Entrez	29817
Expiration Date	Please enquire.