

Product Datasheet

Mouse INSL3 ELISA Kit (Ready to Use) (orb553558)

Description	Mouse Insulin Like Protein 3 (INSL3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	009107
Sensitivity	6.3pg/mL
Entrez	16336
Expiration Date	Please enquire.