

Product Datasheet

Mouse IPF ELISA Kit (Ready to Use) (orb553561)

Description	Mouse Insulin Promoter Factor 1 (IPF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P52946
Sensitivity	0.053ng/mL
Entrez	18609
Expiration Date	Please enquire.