

Product Datasheet

Mouse IL11Ra ELISA Kit (Ready to Use) (orb553588)

Description	Mouse Interleukin 11 Receptor Alpha (IL11Ra) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q64385
Sensitivity	0.061ng/mL
Entrez	16157
Expiration Date	Please enquire.