

Product Datasheet

Mouse IL13Ra2 ELISA Kit (Ready to Use) (orb553592)

Description	Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O88786
Sensitivity	0.126ng/mL
Entrez	16165
Expiration Date	Please enquire.