

Product Datasheet

Mouse KIBRA ELISA Kit (Ready to Use) (orb553641)

Description	Mouse Kidney And Brain Protein (KIBRA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q5SXA9
Sensitivity	30pg/mL
Entrez	211652
Expiration Date	Please enquire.