

Product Datasheet

Mouse KLF1 ELISA Kit (Ready to Use) (orb553650)

Description	Mouse Kruppel Like Factor 1, Erythroid (KLF1) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P46099
Sensitivity	5.8pg/mL
Entrez	16596
Expiration Date	Please enquire.