

Product Datasheet

Mouse LPIN1 ELISA Kit (Ready to Use) (orb553684)

Description	Mouse Lipin 1 (LPIN1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q91ZP3
Sensitivity	0.051ng/mL
Entrez	14245
Expiration Date	Please enquire.