



Product Datasheet

IL17 Receptor B antibody (orb5537)


Description	IL17 Receptor B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human IL-17RB (221-320/502 aa)
Target	IL17RB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q9NRM6
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of rat pancreas tissue using IL17 ...



Immunohistochemical analysis of formalin...



Western blot analysis of rat kidney lysa...