

Product Datasheet

Mouse LPCAT4 ELISA Kit (Ready to Use) (orb553701)

Description	Mouse Lysophosphatidylcholine Acyltransferase 4 (LPCAT4) Sandwich ELISA KitReady to Use: Pre-diluted...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6NVG1
Sensitivity	0.212ng/mL
Entrez	99010
Expiration Date	Please enquire.